

## CARBON AND SULFUR ANALYZER CODE CSA-R300

CAN BE CUSTOMIZED FOR HIGH SULPHUR CHANNEL ACCORDING TO REQUIREMENTS



- Widely used in iron and steel metallurgy, powder metallurgy, magnetic materials, new energy, third-party testing, research institutes, etc.
- Analysis of carbon and sulfur in steel, cast iron, alloys, cement, sand, glass, lime, rubber, catalysts, soil, semiconductors, electronic materials, metal ores, ceramics and other solid materials
- Adoption of high-end infrared detection system and the core components of the air circuit system ensures the stability and accuracy of detection
- The high-frequency heating system adopts special high-frequency capacitors to ensure the long-term reliability of high-frequency work
- Adopt German imported special external metal dust filter, easy to dismantle, easy to clean
- Electronic balance for accurate weighing, auto-matic tare of crucible mass, auto-matic input of sample weight
- The software can monitor solenoid status, signal acquisition status, high frequency status, baseline status, etc.
- With data management and storage, query and statistics, printing projects, generating reports

## STANDARD DELIVERY

Main unit	1pc
Computer	1pc
Software	1pc
Dust catcher	1pc
Electronic balance	1pc
Printer	1pc
Pressure valve	1pc
Standard sample	3pcs
Soda lime (CSA-R300-SL)	1pc
Tungsten flux (CSA-R300-TF)	1pc
Desiccant (CSA-R300-DE)	1pc
Crucible (CSA-R300-CR)	1000pcs
Toolkit*	1set

<sup>\*</sup>Including allen wrench, crucible tweezer, tungsten grain spoon, sample spoon, cotton wool and vacuum silicone grease, etc.

## SPECIFICATION

Analysis accuracy

Sensitivity

Carrier gas

Motive gas

Power supply

Analysis time

Sample weighing

Combustion furnace

Dimensions (L×W×H)

Detector	pyroelectric solid state infrared detectors	
Light source	infrared light source	
Channel configuration	standard with high carbon, low carbon, low sulfur	
Thermostatic	gas chamber thermostat control	
Analysis range	low carbon: 0.6ppm~0.1%, high carbon: 0.1%~50%, low sulfur: 0.6ppm~0.3%*	

recommended range 0.1g~0.5g, can be changed according to the sample content

carbon: 1ppm or RSD  $\leq$  0.5%, sulfur: 1ppm or RSD  $\leq$  1.0%

weight \*The analytical range can be extended to 100% by changing the weighing volume, can be customized for high sulfur channel with 0.3-30% analysis range

oxygen concentration ≥99.5%

nitrogen concentration ≥99.5%

0.01ppm

30S~40S

18MHz, 2.7KVA

AC 220V. 50Hz

70ka

545×380×435mm